

GLUCOSE

Enzymatic-spectrophotometric
GLUCOSE OXIDASE/PEROXIDASE

The application parameters comprised here constitute a guide to facilitate the validation of our reagents by the instrument. It is advisable to validate the use when there is any change in software or reagent versions.

Instrument: CS-T240

Reagent preparation

Reagent is ready to be used.

Instrument settings

Analyze parameters										
Test item	<input type="text" value="Glu"/>	Test full Name	<input type="text" value="Glu"/>	Decimal digit	<input type="text" value="1"/>	Unit	<input type="text" value="mg/dL"/>	<input type="text"/>	Sample blank	<input type="checkbox"/>
Assay	<input type="text" value="1point"/>	Test time	<input type="text" value="10"/>	Point	<input type="text" value="17"/>	<input type="text" value="0"/>	<input type="text" value="0"/>	<input type="text" value="0"/>	Control Interval	<input type="checkbox"/>
Main Wave	<input type="text" value="505"/>	Second Wave	<input type="text" value="0"/>	Instrument factor(Y=aX+b)	a=	<input type="text" value="1.0"/>	b	<input type="text" value="0"/>	Always dilution	<input type="checkbox"/>

Sample Vol.				Reagent					
	Serum			Urine			Vol	Dil	Pos
Normal	<input type="text" value="3"/>	<input type="text"/>	<input type="text"/>	<input type="text" value="3"/>	<input type="text"/>	<input type="text"/>	R1	<input type="text" value="300"/>	<input type="text" value="*"/>
Decrement	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	R2	<input type="text"/>	<input type="text"/>
Increment	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>			
Abs Limit	<input type="text" value="1"/>			<input type="text"/>					Positive reaction
Prozone Limit	<input type="text" value="-3.3"/>			<input type="text"/>					Lower limit

Range parameters

Test item : Glu

Serum		
Age	Male	Female
<input type="text"/>	<input type="text"/>	<input type="text"/>
<input type="text"/>	<input type="text"/>	<input type="text"/>
<input type="text"/>	<input type="text"/>	<input type="text"/>
Years		
Years		
Years		
Default value		
	<input type="text" value="70"/>	<input type="text" value="105"/>
Linear range		
	<input type="text" value="0-500"/>	

Urine		
Age	Male	Female
<input type="text"/>	<input type="text"/>	<input type="text"/>
<input type="text"/>	<input type="text"/>	<input type="text"/>
<input type="text"/>	<input type="text"/>	<input type="text"/>
Years		
Years		
Years		
Default value		
	<input type="text"/>	<input type="text"/>
Linear range		
	<input type="text"/>	

Calibration parametes

Test Item : Glu

		Drift checkup: 3.3
Calibration methods	<input type="text" value="2 point Linearity"/>	Point : <input type="text" value="2"/> Span <input type="text"/>
		Discreteness checkup: 3.3
		Sensitivity checkup: 0
Calibration Solution	Blank horizontal checkup: -3.3-3.3	
(1) 0		
(2) * S1		
* data entered by user		