

LIPASE

Color

The application parameters comprised here constitute a guide to facilitate the validation of our reagents by the instrument. It is advisable to validate the use when there is any change in software or reagent versions.

Instrument: CS-T240

Reagent preparation

Reagent 1: use the Reagent A.
Reagent 2: use the Reagent B.

Analyze parameters										
Test item	<input type="text" value="Lipase"/>	Test full Name	<input type="text" value="Lip"/>	Decimal digit	<input type="text" value="0"/>	Unit	<input type="text" value="U/L"/>	<input type="text"/>	Sample blank	<input type="text"/>
Assay	<input type="text" value="Rate A"/>	Test time	<input type="text" value="10"/>	Point	<input type="text" value="20"/>	<input type="text" value="30"/>	<input type="text" value="0"/>	<input type="text" value="0"/>	Control Interval	<input type="text"/>
Main Wave	<input type="text" value="570"/>	Second Wave	<input type="text" value="0"/>	Instrument factor(Y=aX+b)	a=	<input type="text" value="1.0"/>	b	<input type="text" value="0"/>	<input type="text"/>	Always dilution

Sample Vol.				Reagent			
	Serum		Urine		Vol	Dil	Pos
Normal	<input type="text" value="3"/>	<input type="text"/>	<input type="text" value="3"/>	<input type="text"/>	R1	<input type="text" value="300"/>	<input type="text" value="*"/>
Decrement	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	R2	<input type="text" value="60"/>	<input type="text"/>
Increment	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>			
Abs Limit	<input type="text" value="1"/>		<input type="text"/>				Positive Reaction
Prozone Limit	<input type="text" value="-3.3"/>		<input type="text"/>				Lower limit

Range parameters

Test item : Lipase

Serum									
Age			Male			Female			
Adult		Years							
		Years							
		Years							
Default value									
			0		38				
Linear range									
			5-250						

Urine									
Age			Male			Female			
		Years							
		Years							
		Years							
Default value									
Linear range									

Calibration parameters

Test Item : Lipase

				Drift checkup: 3.3	
Calibration methods	2 point Linearity	Point : 2	Span	Discreteness checkup: 3.3	
Sensitivity checkup: 0					
Calibration Solution				Blank horizontal checkup:-3.3-3.3	
(1)	0				
(2)	* S1				
* data entered by user					