

ALBUMIN (URINE)

Turbidimetry
LATEX

The application parameters comprised here constitute a guide to facilitate the validation of our reagents by the instrument. It is advisable to validate the use when there is any change in software or reagent versions.

Instrument: CS-600 B

Reagent preparation

Reagent 1: Use the Diluent.

Reagent 2: Use the Latex

Instrument settings

Analyze parameters									
Test item	Mau	Test full Name	Mau	Decimal digit	1	Unit	mg/L		Sample blank
Assay	2point end	Test time	10	Point	15	34			Control Interval
Main Wave	546	Second Wave	0	Instrument factor(Y=aX+b)	a=	1.0	b	0	Always dilution

Sample Vol.				Reagent			
	Serum		Urine		Vol	Dil	Pos
Normal	3		3		280		*
Decrement							
Increment							
Abs Limit	1						Positive reaction
Prozone Limit	-3.3						Lower limit

Range parameters

Test item : Mau

Serum			Age		Male		Female	
Adults		Years						
Children		Years						
		Years						
Default value								
			0		15			
Linear range								
			0.9-130					

Urine			Age		Male		Female	
		Years						
		Years						
		Years						
Default value								
Linear range								

Calibration parameters

Test Item : Mau			Drift checkup: 3.3	
Calibration methods	2 point Linearity	Point : 2	Span	Discreteness checkup: 3.3
Calibration Solution			Sensitivity checkup: 0	
(1)	0		Blank horizontal checkup:-3.3-3.3	
(2)	*	S1		

* data entered by user